


Thrive Hearing Control App User Guide

Congratulations on your
journey to better hearing

thrive



Thrive Hearing Control app

Table of Contents

Introduction	4
Apple Users	
Pairing/Connecting	5
Downloading	8
Android Users	
Downloading	9
Pairing/Connecting	10
Home Screen	
Settings	13
Help	13
Connection Status	13
Thrive Assistant	13
Current Memory	13
Thrive Score	13
Volume	14
Mute	14
Navigation	14
Home	15
Memory Menu	15
Customize Screens	
Customize Screens	16
Equalizer	16
Speech in Noise	16
Machine Noise	17
Wind	17
Microphone Direction	18
SoundSpace	18

Table of Contents

Thrive Score Screens	
Overall	19
Brain	19
Body	20
Fall Detection and Alerts	21
Accessories	23
IntelliVoice	23
Settings	
User Guide	24
Device Settings	24
Edge Mode	24
Thrive Score	25
Fall Alert	26
Translate	27
Transcribe	27
Self Check	28
Reminders	28
Hearing Care Anywhere	30
Thrive Care app	31
Find My Hearing Aids	32
App Settings	32
About	32
Demo Mode	33

Features vary by product. This app may have slight differences depending on your phone.

Thrive Compatible Hearing Aids



Shown above are the Thrive platform 2.4 GHz hearing aids that work with the Thrive Hearing Control app.

Apple or Android

The iOS or Android operating system is responsible for managing the connectivity between your Apple or Android device(s) and Thrive compatible hearing aids. The Thrive Hearing Control app is a supplement to your hearing aids offering a variety of enhanced adjustment capabilities. With Thrive, you can customize your hearing aids to fit your lifestyle.

iOS UPDATE NOTE: When updating the version of iOS on your iPhone the hearing aids may lose connection, requiring you to unpair and re-pair the hearing aids with your iPhone.



For the most up-to-date compatibility information, visit starkey.com/thrivesupport

Connectivity for Apple Devices



In order to use the Thrive app, you must first pair your hearing aids with your iOS device.

1. Ensure Bluetooth® setting is enabled on your iOS device. Within the Settings menu go to Bluetooth and toggle to **On**.
2. Turn your hearing aids off and back on. This puts the hearing aids in pairing mode.
3. Within the Settings menu go to **Accessibility > Hearing Devices**.

- You will see your hearing aid name (e.g. "Chris Hearing Aids") when the iOS device first discovers your hearing aids.
- If your name does not appear in the "Devices" list within 5–7 seconds, tap **Accessibility** in the upper left corner, then tap **Hearing Devices**.



4. Tap on the hearing aid name to connect your hearing aids to the iOS device.
5. You will see two pairing requests (one for each hearing aid). Select **Pair** for each request. These requests may be several seconds apart.
6. When pairing is complete, your hearing aid name will change from black to blue.



Reconnecting your hearing aids to your Apple Device

When you turn off your hearing aids or Apple device, they will no longer be connected to one another. To reconnect, turn on your Apple device and then power on your hearing aids.

The hearing aids will automatically connect to your Apple device. You can also verify this connection by accessing the Accessibility shortcut (triple click the Home button) and making sure your hearing aids have been detected.

Download the Thrive app for Apple Devices

The Thrive app makes it easy to control and adjust your hearing aids with your compatible Apple device. Downloading the app is easy. Follow the step-by-step instructions below:

Step One: On a compatible Apple device, open the "App Store."



Step Two: Type "Thrive Hearing Control" in the search box.



Step Three: Tap "GET" to download the app to your device. Follow the screen prompts to set up the app.



Step Four: Once the app has installed, open the Thrive app from your Home screen or Apps List.



Download the Thrive app for Android Devices

The Thrive app makes it easy to control and adjust your hearing aids with your compatible Android device. Downloading the app is easy. Follow the step-by-step instructions below:

Step One: On a compatible Android device, open the "Google Play Store."



Step Two: Type "Thrive Hearing Control" in the search box.



Step Three: Tap "INSTALL" to download the app to your device. Follow the screen prompts to set up the app.



Step Four: Once the app has installed, open the Thrive app from your Home screen or Apps List.



Connectivity for Android Devices

Follow the instructions below to connect your hearing aids to your Android device. If your hearing aids lose connection to your device, repeat the instructions below to reconnect.

1. Launch the Thrive app.
2. After searching, available hearing aids will display on the screen.
NOTE: Each hearing aid will appear separately.
3. Check the box next to each hearing aid and select **Connect**.

NOTE: Pairing can still be completed through the Bluetooth menu of your Android device, if preferred.


Reconnecting your hearing aids to your Android Device

When you turn off your hearing aids or Android device, they will no longer be connected to one another. To reconnect, turn on your Android device and then power on your hearing aids.

The hearing aids will automatically connect to your Android device. If you are having trouble getting your hearing aids and app to reconnect, check the Bluetooth settings on your Android device to make sure your hearing aids are still paired and if not, repeat the pairing process.

Launch/Setup

Upon first launch of the Thrive app, you will be prompted to read and acknowledge the license agreement and privacy policy, as well as customize the setup of the app. You will have the option to set up notification alerts and to set up a cloud account for the Hearing Care Anywhere feature. The Thrive app will default to the Advanced mode; you can choose to set up the Basic mode after installation via the Settings menu.



This app may have slight differences depending on your phone.

Settings

Tap to open the Menu. See Settings section for full description.

Help

The Help feature explains the function of each available screen. All help files can be found in [Settings > User Guide](#).

Connection Status

Tap to view the current connection status for your hearing aids.

Thrive Assistant*

With Thrive Assistant, you can trigger a voice assistant within the Thrive app to ask a question like, "How do I customize my hearing aids settings?" Thrive Assistant uses the phone microphone to listen for your question and then provides you with an answer via a help file or cloud-based services. Thrive Assistant can be triggered via an icon in the Thrive app or a double tap to your Livio Edge AI or Livio premium hearing aids.

Find My Phone helps you find your lost or misplaced smartphone with just a double tap to your Livio Edge AI or Livio premium hearing aids. Use phrases like "I can't find my phone" or "Where is my phone?"

Current Memory

The current memory you're using is shown here.

Thrive Score

A view into both brain and body wellness.

*Requires a cloud account.

Volume


Increase/decrease the volume level on your hearing aids by sliding the volume slider for each side or by using the slidepad (between left and right volume sliders) to adjust both sides simultaneously. Double-tap anywhere on the slidepad to return your volume to the default level set by your clinician.

Mute

Tap the Mute button to mute the microphones on both hearing aids. Tap on it again to unmute the hearing aids.

Navigation

The Navigation bar at the bottom of the screen provides access to the commonly used sections Customize, Thrive Score, Accessories and Settings.




Home

Tap the Home icon on the navigation bar to return to the main screen.

Memory Menu

Tap the Memory Menu icon to view and select from available memories or to create a new memory or edit a memory.




Customize Screens

Tap Customize on the navigation bar to gain access to a variety of adjustment controls for personalizing your hearing aid memories.


Equalizer

Customize offers equalizer bars for several frequency ranges. You can slide up/down to adjust for your preferred sound quality. You can adjust differently per memory.




Speech in Noise

Speech in Noise reduces the amount of background noise coming into your hearing aids. You can adjust this individually for each hearing aid memory. Setting further to the right provides more comfort in noise.




Machine Noise

Machine Noise reduces sounds produced by machinery to increase your comfort in loud environments. Setting further to the right provides more comfort in noise.




Wind

Wind provides you control for reducing the sound of wind noise blowing across your hearing aid microphones per memory. Setting further to the right provides more comfort in wind.




Microphone Direction

Microphone Direction provides the ability to control the sound around you in each memory. Use this option to improve your ability to hear conversation by focusing the microphone to the sounds in front of you. The automatic mode follows speech sounds based on noise levels.



SoundSpace

SoundSpace allows you to turn your iOS device into a personal settings adjustment tool for your hearing aids. It provides the ability to modify the response for any of your custom Thrive memories to better suit a particular environment. As you move your finger around the blank screen, the settings in your hearing aids will temporarily change so you can listen to a variety of options.



Thrive Score


Overall Thrive Wellness Score is a combination of both Brain and Body scores. You have the ability to reach a score of 200 each day.

Brain tracking is comprised of three segments. Each segment has a point value to reflect how you are using your hearing aids. The Thrive app reads usage data from your hearing aids and shows your daily progress points. Your goal is 100 points per day.

Use Score points are gained by wearing your hearing aids on a daily basis. The more you use them, the more points earned (40 possible).

Engagement Score points are gained by using your hearing aids to engage with others. The more time spent in environments with conversational speech, the more points earned (40 possible).

Environment points are gained when the hearing aids are adapting to challenging listening environments (20 possible).




Body* tracking is comprised of progress toward three daily goals for physical activity. Meeting all three goals will provide a full 100 points. The hearing aid tracks steps and that data is transmitted to the app throughout the day. You can set/modify goals in the Settings area.

Steps shows your progress toward the daily step count goal (40 points possible).

Exercise shows your progress toward your daily goal for vigorous exercise (40 points possible).

Stand shows your progress toward your move goal to stand up and move around every hour (20 points possible).




NOTE: Hearing aids must regularly connect to the phone during the day and the Thrive app must be kept open in the background of your phone.

**For Livio Edge AI hearing aids.*

Fall Detection and Alerts*

The Fall Alert feature can be used to send an alert text message containing your GPS location to your designated contacts, should you fall or experience a non-fall-related event. Fall Alert can be configured to send automated and/or manually initiated alerts.




Auto Alert enables the sensors in your hearing aids to detect a fall automatically.

Manual Alert enables you to “Push and Hold” the user control on either hearing aid to initiate an alert text message.

An alert text message will be initiated by the Thrive app when a fall is automatically detected or manually initiated by you. The text message will be sent to a maximum of three predefined contacts, who can confirm receipt of the message and view a map, indicating your location.

**For Livio Edge AI hearing aids.*




Alert Cancellation allows you to cancel an Auto Alert or Manual Alert from either your hearing aids or your smartphone. To cancel an alert text message from your hearing aids, press the user control on either hearing aid. You will be informed of the cancellation in the Thrive app and via a speech indicator through your hearing aids. Fall Alert messages can be canceled within the 60 or 90 second preselected cancellation time following alert initiation.


Contacts allows you to identify up to three contacts to whom you would like alert text messages sent. You must enter the name and smartphone number for each contact into the Thrive app. Each of your contacts will receive a text message prompting them to confirm participation in your Fall Alert system.

Accessories

The Thrive app provides the ability to control streaming audio from most 2.4 GHz accessories. These accessories pair directly to your hearing aids. By tapping on the Accessories icon within the Navigation bar, you can select a paired accessory, start and stop audio streaming and control the volume of both streaming and surrounding sounds.



IntelliVoice* combines the advantages of a microphone in your smartphone with the computational power of artificial intelligence. The benefit is more powerful and flexible speech enhancement in the most demanding listening environments for patients with severe hearing loss only. Ask your clinician for more information.



**For Livio Edge AI hearing aids.
Available with iOS smartphones only.*

User Guide

An overview of each screen/feature is available in the Thrive app. Each section shows the screen details and describes each function.

Device Settings

Advanced settings for your hearing aids can be turned on/off or set in this area.

Tap Gesture allows you to enable or disable the tap gesture, set the sensitivity, and choose the feature to activate by tap gesture (Thrive Assistant, Accessory Start/Stop Streaming or Edge Mode**).

Thrive Assistant provides you the ability to use your voice to ask questions about Thrive, ask general questions that are answered from the Internet (if you've logged in), command your hearing aids to make changes to memory or volume, or set a reminder. Examples for voice command: "change my memory to restaurant" or "set a reminder to pick up my prescription tomorrow at 3:00."

Edge Mode* allows you to trigger an automatic assessment of the current listening environment and Edge AI hearing aids will create optimized settings for you based on the characteristics of that particular environment. These settings are in place until you change memories or restart your hearing aids.

**For Livio Edge AI hearing aids.*

Auto Streaming enables your hearing aids to recognize when an audio stream starts and automatically change to a different memory and microphone reduction specifically optimized for streaming audio.

Comfort Boost allows your hearing aids to aggressively reduce noise while optimizing sound quality in noisy environments.

Alerts allows your hearing aid to play an alert tone for certain phone notifications like calendar events, email and social media.

Auto On/Off easily and automatically saves battery power when you're not using your hearing aids. Take them off, set them on a flat, stable surface, and they will go into a low-power mode after 15 minutes. Pick your hearing aids up and put them on again to automatically resume full power.

Thrive Score

Settings section for the Thrive Score. You can add/update your demographic data for more accurate step-tracking activity. In addition, you can set/modify goals for Steps, Exercise and Stand.

Fall Alert*

Auto Alert sensitivity allows you to adjust your Auto Alert sensitivity in the Thrive app. Increasing the sensitivity may increase the likelihood of detecting a fall. Decreasing the sensitivity may help reduce the probability of false alerts.

Profile enables you to modify your profile information.

Auto Alert and Manual Alert indicators provide the ability to review and demo the Fall Alert indicators for familiarity in the event of a fall.

Contacts enables you to add, edit or delete a contact.

Cancellation Timer allows you to select 60 or 90 seconds as the allotted time before an alert text message is sent.

Alert Message enables you to edit your name and displays an example of an alert text message.

*For Livio Edge AI hearing aids.



Translate*

Translate is a feature that allows you to choose your language and the language of another speaker to engage in conversation. When you speak into the phone, the app will translate your speech and display in the other person's language. When the other person speaks into the phone, the app will translate their speech, display it in your language and stream that translated text to your hearing aids in your language.

Transcribe*

Transcribe is a feature that allows you to have closed captioning for real world situations. When you or others speak into the phone, the Thrive app will transcribe that speech and display it onto your screen. There is also the option to save, message, copy or email the transcribed text.

*For Livio Edge AI hearing aids.





Self Check

Self Check provides you a quick, convenient way to analyze your hearing aids' performance on your own. With Self Check, you are able to run diagnostics on the performance of the components within your hearing aid(s) including the microphones, receiver, circuit and sensors.

Self Check will display the performance of components and provide you with instructions for servicing if needed (most issues are due to wax/debris and can be resolved easily). For AI and regular hearing aids, your hearing professional will do the initial setup of Self Check to establish the baseline measure for you. For Edge AI hearing aids, the baseline will come already established.

Reminders



The Reminders feature can help you keep track of important tasks and appointments in your busy life, without having to remember them all.



You can set up and receive visual reminder notifications on your smartphone and/or audible reminder notifications streamed from your smartphone through your hearing aids for a multitude of tasks, such as taking medications, refilling prescriptions, cleaning your hearing aids, exercising and drinking water.

In addition, you can customize your reminders, i.e., Take blood pressure medication, Walk the dog, Call the physical therapist, Bring water bottle and mat to yoga. After you set the reminder name, date, and time, you can select the option to have this reminder repeat daily, weekly, monthly or yearly.

And you need not worry if you are unable to complete the task the moment you receive your reminder, because overdue tasks will be displayed in red on the Reminders screen until you check them off.



Hearing Care Anywhere

Creating an Account


Hearing Care Anywhere requires you to create your own account. Through this account, your hearing aids can be backed up and synched between multiple phones. You can create an account during the install process or at a later time if desired. You can create an account using your existing Facebook or Google account, or by creating your own account using email and a password of your choosing.

Back Up

When you make any changes using the Thrive app, those changes are saved to your cloud account to back up your settings in real-time.

Live Sessions


Your clinician can enable Live Sessions. This allows your clinician to use video chat for hearing aid adjustments in real-time without having to stop into the office.



Thrive Care app*

The Thrive Care app provides peace of mind to those you choose to share helpful information with while helping you feel confident and independent.

With Thrive Care, you can share information like physical activity, hearing aid usage, social engagement and more, with your caregivers or loved ones.



*For Livio Edge AI hearing aids.

Find My Hearing Aids


The Find My Hearing Aids feature allows you to track the last known physical location of your hearing aids on a map. When the Find My Hearing Aids feature is accessed, you will see a map view with a note showing the location where your hearing aids were last connected. If your hearing aids are actively connected (or become connected), the bars on the bottom of the screen will animate to indicate how close or far away you are from your hearing aids.

App Settings

Options under App Settings include the ability to change between Basic or Advanced mode. The Advanced mode provides full functionality for every feature in the app, while the Basic mode includes a subset of just the very basics to keep things simple.

About

Provides information about Starkey, the End User License Agreement, privacy policy and the app version. There is also a section to provide feedback.



Demo Mode

Puts the Thrive app into a demo mode allowing the app to function without any hearing aids attached so that all features are visible for demonstration. Adjustments made while in Demo Mode will not be applied to your hearing aids.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

For Thrive Hearing Control app and iOS/Android compatibility information, please visit starkey.com/thrivesupport.

Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a registered service mark of Apple Inc.

Android, Google and Google Play are trademarks of Google LLC.

The Starkey logo, Livio, Hearing Care Anywhere, SoundSpace, Thrive and Thrive logo are registered trademarks of Starkey Laboratories, Inc. Thrive Care, Thrive Care logo and IntelliVoice are trademarks of Starkey Laboratories, Inc.

iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Starkey is under license.

The
**American
Hearing
Company**

Veteran technical support line*:

(800) 211-7789

*For technical support of hearing devices with wireless connectivity



Starkey.com/veterans

6700 Washington Ave. S.
Eden Prairie, MN 55344
1.866.418.6928

