


Thrive 2.4 GHz + NFMI* TECHNOLOGY

livio Edge AI

livio AI

livio



2.4 GHz ACCESSORIES



CROS Transmitter

Our wireless CROS System takes audio from an ear-level microphone and wirelessly transmits it to the opposite ear. Available for RIC 312, RIC Rechargeable and BTE styles.



Table Microphone

Enjoy group settings like family gatherings or noisy restaurants with eight built-in microphones that detect the primary speaker's voice and stream it to 2.4 GHz hearing aids.



Remote Microphone +

Enjoy one-on-one conversations in noisy environments and easily stream binaurally from four different audio sources: Bluetooth®, Telecoil, FM or Line-In.



Mini Remote Microphone

Patients with older mobile devices can enjoy one-on-one conversations in noisy environments. Clip it onto the clothing of the person being spoken to or use it as a TV streamer by placing it near the sound source.



Remote

Also ideal for patients using older mobile devices with updated features so they can control memory and volume, mute their 2.4 GHz hearing aids and turn other special features on and off.



Standard & Custom Charger

Both chargers hold enough power to provide portable charging three times. LED indicators let patients know when their hearing aids are fully charged, and hearing aids turn on automatically when removed from the charger.



Mini Turbo Charger

A pocket-sized lithium-ion solution that provides a 3.5-hour charge in just seven minutes. A fully charged Mini Turbo Charger can completely charge a pair of RIC R or BTE R hearing aids four times.

Livio Edge AI Premium 2400	Livio AI Advanced 2000	Livio AI Select 1600	Livio AI Basic 1200
----------------------------	------------------------	----------------------	---------------------

ARTIFICIAL INTELLIGENCE FEATURES	Healthable Technology:			
	Fall Alert	●		
Body Tracking	●	●	●	●
Brain Tracking	●	●	●	●
Intelligent Assistant:				
Translate	●			
Transcribe	●	●		
Thrive Assistant	●	●	●	
Reminders	●	●	●	
Find My Phone	●	●	●	
Auto On/Off	●	●	●	●
Self Check	●	●	●	●
Tap Control	●	●	●	●

Platform:	Thrive	Thrive	Thrive	Thrive	Thrive	Thrive
Sound Imaging: Channels Bands	24	24	20	16	12	10
Edge Mode	●					
IntelliVoice	●					
Speech Optimization	●	●	●	●	●	●

HEARING REALITY PRO*	Music Optimization:					
	Music Adaptation	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
E2E Music Adaptation	●	●	●	●	●	●

HEARING REALITY PRO*	Ear-to-Ear Technology:					
	E2E Wind Noise Management	●	●	●	●	
E2E Machine Noise Adaptation	●	●	●	●		
E2E Directionality***	●	●	●	●		
E2E Phone Streaming*	●	●	●	●		

HEARING REALITY PRO*	Sound Manager:					
	Auto Music	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Quiet	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Spatial Speech Enhancement	●	●	●	●		
Machine Noise	Up to 20 dB of reduction	Up to 20 dB of reduction	Up to 10 dB of reduction	Up to 5 dB of reduction	Up to 5 dB of reduction	Up to 5 dB of reduction
Speech in Noise	Up to 20 dB of reduction	Up to 20 dB of reduction	Up to 8 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction
Transient Noise Reduction	Up to 15 dB of reduction	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction
Wind	Up to 30 dB of reduction	Up to 30 dB of reduction	Up to 15 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction

Livio Edge AI Premium 2400	Livio Premium 2400	Livio AI/Livio Advanced 2000	Livio AI/Livio Select 1600	Livio AI/Livio Basic 1200	Livio Economy 1000
----------------------------	--------------------	------------------------------	----------------------------	---------------------------	--------------------

HEARING REALITY PRO*	Directional Processing:					
	Immersion	●	●	●	●	●
Adaptive	●	●	●	●	●	
Dynamic	●	●	●	●	●	●
Directional	●	●	●	●	●	●
Feedback Management	●	●	●	●	●	●
Frequency Lowering	●	●	●	●	●	●
Tinnitus Technology	●	●	●	●	●	●
CROS System	●	●	●	●	●	
Telecoil	●	●	●	●	●	●
2.4 GHz Accessory Compatibility	●	●	●	●	●	●

HEARING REALITY PRO*	Telehealth Service:					
	Hearing Care Anywhere (Live Sessions & Remote Adjustments)	●	●	●	●	●



thrive & thriveCare

Give patients more control! Fully updated for improved usability and overall performance – the **Thrive Hearing Control app** is better than ever. Peace of mind! The **Thrive Care app** offers Livio Edge AI patients reassurance and independence.

*Available on RIC R, RIC 312 and BTE 13 styles only. | **Hearing Reality Pro is exclusive to Livio Edge AI. | ***Available for products with directional microphones only.

Synergy 900SYNC TECHNOLOGY

Muse iQ



Power Plus BTE 13

P **T** **C** **TC**

Picasso



ITE

T **TC**



ITC/HS

T **TC**



CIC

T



IIC*

T

Feature	Premium i2400	Advanced i2000	Select i1600	Basic i1200	Economy i1000
Platform:	Synergy	Synergy	Synergy	Synergy	Synergy
Sound Imaging: Channels Bands	24	20	16	12	10
Acuity Speech Optimization:	●	●	●		
Music Optimization:					
Music Adaptation	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
E2E Music Adaptation	●	●	●	●	●

Feature	Premium i2400	Advanced i2000	Select i1600	Basic i1200	Economy i1000
Ear-to-Ear Technology: Acuity Binaural Imaging					
E2E Wind Noise Management	●	●	●		
E2E Machine Noise Adaptation	●	●	●		
E2E Directionality**	●	●	●		
E2E Phone Streaming	●	●	●		
Environment Manager: Acuity Lifescape Analyzer					
Auto Music	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Quiet	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Machine Noise	Up to 20 dB of reduction	Up to 10 dB of reduction	Up to 5 dB of reduction	Up to 5 dB of reduction	Up to 5 dB of reduction
Speech in Noise	Up to 20 dB of reduction	Up to 8 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction	Up to 6 dB of reduction
Wind	Up to 30 dB of reduction	Up to 15 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction

Feature	Premium i2400	Advanced i2000	Select i1600	Basic i1200	Economy i1000
Directional Processing*:					
Omni	●	●	●	●	●
Adaptive	●	●	●	●	
Dynamic	●	●	●	●	●
Directional	●	●	●	●	●
Feedback Management	●	●	●	●	●
Frequency Lowering	●	●	●	●	●
Tinnitus Technology	●	●	●	●	●
CROS System	●	●	●	●	●
Telecoil	●	●	●	●	●
SurfLink Accessory Compatibility	●	●	●	●	●

P Absolute Power/Power Products

T Tinnitus Technology

C CROS System

TC Telecoil

SURFLINK ACCESSORIES



SurfLink Remote

Patients can control memory and volume adjustments and even adjust their hearing aids, or easily go in and out of streaming mode.



SurfLink Media 2

A set-and-forget media streaming solution that connects to TVs and MP3 players to wirelessly stream audio directly to wireless Picasso and Muse iQ hearing aids.



SurfLink Remote Microphone 2

Stream audio directly to Muse iQ, and wireless Picasso hearing aids or use to help patients in one-on-one conversations or settings with multiple speakers.



SurfLink Mini Mobile System

Patients can connect any device that's compatible with Bluetooth technology to wireless Muse iQ or Picasso hearing aids.

*IIC available in i2400/i2000/i1600 technology tiers only. | **Available for products with directional microphones only. | ***Available for Muse iQ Power Plus BTE 13 style only.

Starkey logo, Thrive, Thrive logo, Livio, Livio logo, Hearing Reality, Hearing Care Anywhere, Thrive Care, Thrive Care logo, IntelliVoice, Muse, Picasso, Synergy and Acuity are registered trademarks of Starkey Laboratories, Inc. Muse logo, and Picasso logo are trademarks of Starkey Laboratories, Inc.

Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Starkey is under license.