

















# Earmolds Reference Guide












## Materials & Options

### Custom RIC Earmolds

Hard Material Styles DIGITAL SLS/LUCITE				Soft Material Styles HIGH STRENGTH SILICONE			
							
<b>RIC EARMOLD</b> (IN LUCITE)*	<b>RIC HELIX</b>	<b>RIC CANAL LOCK</b> (IN LUCITE)*	<b>RIC SKELETON</b> (IN LUCITE)*	<b>RIC HOLLOW</b>	<b>RIC HELIX</b>	<b>RIC EARMOLD</b>	<b>RIC SKELETON</b>
							
<b>AP EMBEDDED</b>	<b>RIC SHELL</b> (IN LUCITE)*	<b>RIC SEMI SKELETON</b> (IN LUCITE)*	<b>RIC HALF SHELL</b> (IN LUCITE)*	<b>RIC CANAL LOCK</b>	<b>RIC SHELL</b>	<b>RIC HALF SHELL</b>	<b>RIC SEMI SKELETON</b>

DEFAULT COLOR FOR BOTH MATERIALS IS CLEAR

### Custom Thin Tube Earmolds

Hard Material Styles DIGITAL SLS/LUCITE				Soft Material Styles HIGH STRENGTH SILICONE			
							
<b>SHELL</b>	<b>HALF SHELL</b>	<b>CANAL</b>	<b>CANAL LOCK</b>	<b>SHELL</b>	<b>HALF SHELL</b>	<b>CANAL</b>	<b>CANAL LOCK</b>
							
<b>SKELETON</b>	<b>HOLLOW CANAL LOCK</b>			<b>SKELETON</b>			











DEFAULT COLOR FOR BOTH MATERIALS IS CLEAR

\*Will default to hollow when ordered in lucite material

# Earmolds Reference Guide

## Materials & Options

### Custom BTE Earmolds

Hard Material Styles DIGITAL SLS/LUCITE				Soft Material Styles HIGH STRENGTH SILICONE			
							
SHELL	HALF SHELL	CANAL	CANAL LOCK	SHELL	HALF SHELL	CANAL	CANAL LOCK
							
SKELETON				SKELETON			

DEFAULT COLOR FOR BOTH MATERIALS IS CLEAR

### Earmold Vent Size Chart

VENT NAME	APPROXIMATE SIZE
Pressure Vent	0.8mm
Small	1.2mm
Medium	1.6mm
Large	2.55mm
Select-A-Vent	up to 3.37mm
Open Fit	As large as canal allows

VENT SIZE IS LIMITED BY EAR ANATOMY