

It's time.

The ultimate hearable meets the ultimate wearable.

Halo™ 2 – compatible with Apple Watch®

Another hearing breakthrough that will help attract more patients to your practice ... Halo 2 hearing aids and Apple Watch come together to create happier, more connected patients.

Halo™ 2



TruLink™ Hearing Control app for Apple Watch provides your patients even more ways to easily control and personalize their hearing experience – so they can hear better and live better.

In addition to all of the TruLink functionality with iPhone®, patients can use their Apple Watch to:

- ▶ **Control volume easily and effortlessly** in any environment
- ▶ **Change memories to personalized settings** to optimize audibility, comfort and sound quality
- ▶ **Mute hearing aid microphones** to block out surrounding noises

Made for
iPod iPhone iPad

Download on the
App Store



Hearing Is Our Concern™



Tinnitus Integration

Integrated controls in the app help patients control the volume and modulation rate of our MultiFlex Tinnitus Technology.



Real-time Notifications

Any feature like new mail, calendar alert, tweet, Facebook post that is set to have an audible alert in Settings/Sound on your phone will now alert your hearing aids to play a unique indicator.



SoundSpace

TruLink SoundSpace enables patients to easily adjust sound settings to their specific environment. By simply moving their finger on the screen, they can adjust frequency and tone to suit their needs. They're also able to save and fine-tune their preferences in different situations as a TruLink Memory so that they're comfortable in any listening environment*.



Remote Control

With the touch of their smartphone, patients can easily change volume and switch Memories on their hearing aids. Both are adjustable with a visual indicator bar for better control and nearly infinite levels to suit their preference.



TruLink Memories

Patients can create up to 20 TruLink Memories using the SoundSpace tool on the TruLink app. Each memory is programmed and run via the app, and gives patients the ability to label them to their liking. Patients can even geotag their TruLink Memories using the built-in GPS on their phone. A geotagged memory will recognize where they are, and automatically adjust their Halo hearing aids when they enter that space. For instance, a geotagged "coffee shop" memory will automatically activate when they walk into their favorite coffee shop.



TruLink Beacons

A TruLink beacon is a small, low-energy Bluetooth® device that broadcasts a unique identifier your Apple® or Android™ device and TruLink app can recognize and use to automatically change memories in the hearing aid. Beacons can be placed in any location where you want an automatic memory change. For example, place a beacon in any room of your house to auto-change to a memory configured for that room [noisy kitchen, watching television].



Find My Hearing Aids

Patients can easily find lost hearing aids using the Find My Hearing Aids feature within the TruLink app. The Find My Hearing Aids feature shows both a location and a time stamp. A signal detector sends a stronger or weaker signal based on how close patients are to locating their hearing aids.



Apple Watch Compatible

Use the TruLink app for Apple Watch to control volume, change memories and mute the patients' hearing aids right from their wrists.

*SoundSpace not available for new Music Memory for Halo 2 due to the complex nature of music processing.

"Made for iPod," "Made for iPhone," and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, iPad or Apple Watch may affect wireless performance. Apple, the Apple logo, iPhone, iPad, iPod touch, Apple Watch and App Store are trademarks of Apple Inc., registered in the U.S. and other countries. Halo 2 and TruLink are compatible with iPhone 6s Plus, iPhone 6s, iPhone 6 Plus, iPhone 6, iPhone 5s, iPhone 5c, iPhone 5, iPhone 4s, iPad Air 2, iPad Air, iPad (4th generation), iPad mini 3, iPad mini with Retina display, iPad mini, iPod touch (5th generation) and Apple Watch.